ZIRKOFAST PLUS

SINTERING FURNACE FOR RAPID ZIRCONIA

The **ZIRKOFAST PLUS** furnace is designed for the sintering of both rapid and standard zirconia for dental prostheses and is equipped with an advanced electronic control system featuring a microprocessor and permanent memory. It allows the storage of up to 50 heating programs, each consisting of up to 6 phases (4 heating phases and 2 cooling phases), with adjustable temperature, heating rate in °C/min, and dwell times.

All cycles are fully customizable and can be modified by the operator at any time. The furnace can perform cycles ranging from a minimum heating rate of 1 °C/min (with a total heating duration of approximately 26 hours) to a maximum heating rate of 60 °C/min to reach the working temperature of 1.650°C in about 30 minutes only for the rapid Zirconia, with the possibility of setting different rates for each phase.

The furnace controls both heating and cooling, allowing dwell phases even during cooling. The cooling rate can also be adjusted, noting that the maximum cooling speed corresponds to natural cooling, while the minimum can be set freely, similar to the heating phases.


The MoSi₂ heating elements ensure reliable operation up to 1800 °C. They are positioned in view inside the chamber and heat by radiation, with a maximum supply voltage of 24 V.

The door is equipped with a mechanical safety lock and a micro-switch that automatically interrupts power to the heating elements when the door is opened and the furnace features a Touch Screen display for simple and intuitive user interaction.

USB Port: Programs can be stored or uploaded via USB flash drive, allowing convenient sharing and transfer of cycles, including via e-mail.

Reliable and safe is Made in Italy under CE normatives with 12 months in-house warranty.

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The support is stackable so it is possible to insert up to 3 stacked pieces into the chamber.

TISSIDental

Dental and Laboratory Equipment Made in Italy

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